


# Piezo Positioners

## Short Instructions



P-5xx / P-6xx / P-7xx



### User information

These short instructions contain an overview of the most important safety and handling instructions for piezo positioners with the product codes given above (x: any number).

Subject to change. These short instructions are superseded by any new release. The latest respective release is available for download on the PI website.



### Downloading and reading the manual

The actions during installation, startup, operation, and maintenance require additional information from the manuals of your product. Manuals can be titled as follows: "User Manual", "Technical Note".

#### Downloading manuals from the website

1. Open the website [www.pi.ws](http://www.pi.ws).
2. Search the website for the product number (e.g. P-725).
3. Select the corresponding product to open the product page.
4. Select the **Downloads** tab.  
→ The manuals are shown under **Documentation** .
5. Select **ADD TO LIST** for the desired manual and then select **REQUEST**.
6. Fill out the request form and select **SEND REQUEST**. The download link will then be sent to the entered e-mail address.

If you cannot find the manual you are looking for or if you have any questions: Contact the PI customer service via [service@pi.de](mailto:service@pi.de).



## Safety

### Intended use

This product is a laboratory device as defined by DIN EN 61010-1. The product is intended for indoor use and use in an environment that is free of dirt, oil, and lubricants.

The positioner is designed for fine positioning as well as moving smaller objects quickly and precisely. The positioner is not intended for applications in areas where failure would be a considerable risk for people or the environment.

It is only possible to use the product as intended when completely installed and connected, and in conjunction with suitable electronics.

The product may only be installed, operated, maintained, and cleaned by authorized and appropriately qualified personnel.

### Electrical dangers

**CAUTION!** Temperature changes and compressive stress can induce charges in the positioner. After disconnection from the electronics, the piezo actuators can remain charged for several hours. Touching charged parts can result in minor injuries from electric shock.

- Do not touch the contacts in the connector.
- Before mounting, connect the positioner to the switched-off electronics for a few seconds to discharge it.
- Secure the connectors with screws against being pulled out.
- Do not pull connectors out of the electronics during operation.

**CAUTION!** If a protective earth conductor is not properly connected, touching the positioner in the case of malfunction can result in minor injuries from electric shock.

- Operate the positioner only with a properly connected protective earth conductor.

- Do not remove the protective earth conductor from the positioner during operation.
- Pay attention to the applicable standards for attaching the protective earth conductor.

**NOTICE!** Moisture, liquids, and electrically conductive materials (e.g., metal dust) that penetrate the positioner can destroy the positioner's piezo actuator.

- Operate the stage only under the permissible ambient conditions (see manual).
- Do not operate a vacuum-compatible positioner during evacuation.

**NOTICE!** Excessively high or wrongly connected operating voltages can cause damage to the positioner.

- Pay attention to the operating voltage range of the positioner (see manual).
- Pay attention to the correct pin assignment.

**NOTICE!** Continuously high voltages can reduce the lifetime of the positioner.

- Keep the maximum operating voltage as low as possible during operation.
- Discharge the positioner appropriately when it is not in use.

### Mechanical dangers

**NOTICE!** Mechanical forces can destroy the positioner.

- Avoid shocks and drops.
- Pay attention to the maximum permissible forces (see manual).
- Do not pull the connecting cable of the positioner out of the electronics during operation.
- If the connecting cable is accidentally pulled out of the positioner during operation: Switch off the electronics before you reconnect the positioner.

- Ship the positioner in the original packaging only.

**NOTICE!** Incorrectly fitted screws can cause damage.

- Do not allow screw heads to protrude.
- Select the screw length according to the depth of the mounting holes.

**NOTICE!** Collisions can damage the positioner, the load to be moved, and the surroundings.

- Mount the positioner and the load so that the load cannot get jammed or blocked, or collide with objects in the workspace.

**NOTICE!** Unsuitable mounting could warp the positioner and reduce the accuracy.

- Mount the positioner on a flat surface with similar thermal expansion properties (see manual for recommended flatness).

**NOTICE!** Uncontrolled oscillation can damage your application or the positioner.

- If oscillation occurs, switch off the servo mode or stop the positioner immediately.
- Check the settings of the servo control parameters (see manual).

## Handling of vacuum-compatible products

**NOTICE!** If the piezo actuator is vacuum-compatible, attention must be paid to cleanliness.

- Only touch the piezo actuator with powder-free gloves.
- If necessary, wipe the piezo actuator clean.




## Installing


### Removing the transport safeguard

If the positioner was delivered with a transport safeguard:

1. Lay the positioner onto a surface with the supporting frame facing downwards.
2. Unscrew the four plastic screws.




3. Remove the plastic screws and the supporting frame.



## Mounting the positioner and affixing the load

1. Mount the positioner to the provided mounting holes.
2. Affix the load to the provided mounting holes.
3. Check that the positioner and the load are affixed firmly.


### Correct mounting



## Mounting on uneven surface



## Incorrect alignment of the load



## Connecting the positioner to the protective earth conductor


### Suitable protective earth conductor

- Cross section of the protective earth conductor  $\geq 0,75 \text{ mm}^2$ .

### Connecting via separate protective earth connector

If a separate protective earth connector is available, it must be used.

1. Attach a suitable cable lug to the protective earth conductor.
2. Use the screw set supplied to attach the cable lug of the protective earth conductor to the protective earth connector.
3. Tighten the M4 screw with at least three turns and a torque of 1.2 to 1.5 Nm.



- 1 M4 threaded hole
- 2 Safety washer
- 3 Flat washer
- 4 Cable lug of the protective earth conductor
- 5 M4 screw

### Connecting via mounting holes

If a separate protective earth connector is not available, the stage has to be mounted on an electrically conductive surface connected to a protective earth conductor.

1. Place a safety washer on each mounting screw.
2. To fix the stage to underlying surface, use sufficiently conductive screws in all mounting holes.
3. Tighten the mounting screws with at least three turns (for torque, see manual).




## Old equipment disposal

In accordance with EU law, electrical and electronic equipment may not be disposed of in EU member states via the municipal residual waste. This refers to devices which are labeled with the symbol for the separate collection of waste electrical and electronic equipment.

If you are established in a member state of the European Union and have an old device from PI, Physik Instrumente (PI) SE & Co. KG will take care of the environmentally friendly disposal of this old device free of charge. You can send the old PI device free of charge to the following address:

Physik Instrumente (PI) SE & Co. KG  
 Auf der Römerstraße 1  
 76227 Karlsruhe  
 Germany

If you are not established in a member state of the European Union: Dispose of your old device according to the national and local rules and regulations.



- 1 Mounting screw
- 2 Safety washer
- 3 Mounting hole



## Startup

The instructions required for startup and operation are in the manual for the electronics used.

## Companies

### Physik Instrumente (PI) SE & Co. KG

Auf der Römerstraße 1  
 76228 Karlsruhe  
 Germany  
 Phone: +49 721 4846-0  
 Fax: +49 721 4846-1019  
 info@pi.de  
 www.pi.ws

### PI miCos GmbH

Freiburger Straße 30  
 79427 Eschbach  
 Germany  
 Phone: +49 7634 5057-0  
 Fax: +49 7634 5057-99  
 info@pimicos.de  
 www.pi.ws

### PI Ceramic GmbH

Lindenstraße  
 07589 Lederhose  
 Germany  
 PI Ceramic GmbH  
 Phone: +49 36604 882-0  
 Fax: +49 36604 882-4109  
 info@piceramic.de  
 www.piceramic.com