

PZ255E
S-334 Piezo Tip/Tilt Mirror
User Manual

Version: 1.0.2

Date: 10.10.2016



This document describes the following products:

- **S-334.2SD1**
High-DynamicsPiezo Tip / Tilt System, 50 mrad,
SGS, Sub-D Connector, incl. Mirror
- **S-334.2SL1**
High-DynamicsPiezo Tip / Tilt System, 50 mrad,
SGS, LEMO Connectors, incl. Mirror



The following company names and brands are registered trademarks of Physik Instrumente (PI) GmbH & Co. KG:

PI®, NanoCube®, PICMA®, PIFOC®, PILine®, NEXLINE®, PiezoWalk®, PicoCube, PiezoMove, PIMikroMove, NEXACT®, Picoactuator®, Plnano®, NEXSHIFT, PITOUCH, PIMag®, PIHera, Q-Motion®

© 2016 Physik Instrumente (PI) GmbH & Co. KG, Karlsruhe, Germany. The text, photographs and drawings in this manual are protected by copyright. With regard thereto, Physik Instrumente (PI) GmbH & Co. KG retains all the rights. Use of said text, photographs and drawings is permitted only in part and only upon citation of the source.

Original instructions

First printing: 10.10.2016

Document number: PZ255E, KSch, Version 1.0.2

Subject to change without notice. This manual is superseded by any new release. The latest release is available for download (p. 3) on our website.

Contents

1	About this Document	1
1.1	Objective and Target Audience of this User Manual.....	1
1.2	Symbols and Typographic Conventions.....	1
1.3	Figures	2
1.4	Other Applicable Documents	2
1.5	Downloading Manuals.....	3
2	Safety	5
2.1	Intended Use	5
2.2	General Safety Instructions	5
2.3	Organizational Measures.....	6
3	Product Description	7
3.1	Model Overview	7
3.2	Product View	8
3.3	Product Labeling.....	10
3.4	Scope of Delivery.....	13
3.5	Accessories	13
3.6	Suitable Electronics	14
3.7	Control.....	15
3.8	Mirror	16
4	Unpacking	17
5	Installation	19
5.1	General Notes on Installation.....	19
5.2	Mounting the S-334 and Connecting it to a Protective Earth Conductor	21
5.3	Removing the Transport Lock.....	23
5.4	Connecting the S-334 to the Controller	24
6	Start-Up and Operation	25
6.1	General Notes on Start-Up and Operation.....	25
6.2	Operating the S-334	27
6.3	Discharging the S-334.....	28

7	Maintenance	29
7.1	General Notes on Maintenance	29
7.2	Cleaning the S-334.....	30
7.3	Packing the S-334 for Transport.....	30
8	Troubleshooting	33
9	Customer Service	35
10	Technical Data	37
10.1	Specifications.....	37
10.1.1	Data Table.....	37
10.1.2	Maximum Ratings.....	38
10.1.3	Ambient Conditions and Classifications	39
10.2	Dimensions	40
10.3	Pin Assignment	41
10.3.1	S-334.2SD1: Sub-D 25 (m) Piezo and Sensor Connection.....	41
10.3.2	S-334.2SL1: LEMO Piezo and Sensor Connections	42
11	Old Equipment Disposal	43
12	EC Declaration of Conformity	45

1 About this Document

In this Chapter

Objective and Target Audience of this User Manual	1
Symbols and Typographic Conventions	1
Figures.....	2
Other Applicable Documents.....	2
Downloading Manuals	3

1.1 Objective and Target Audience of this User Manual

This user manual contains the information needed for the intended use of the S-334.

Basic knowledge of servo systems, drive technologies, and suitable safety measures is assumed.

The latest versions of the user manuals are available for download (p. 3) on our website.

1.2 Symbols and Typographic Conventions

The following symbols and typographic conventions are used in this user manual:

CAUTION



Dangerous situation

If not avoided, the dangerous situation will result in minor injury.

- Actions to take to avoid the situation.

NOTICE




Dangerous situation

If not avoided, the dangerous situation will result in damage to the equipment.

- Actions to take to avoid the situation.

INFORMATION

Information for easier handling, tricks, tips, etc.

Symbol/ Label	Meaning
1. 2. ➤	Action consisting of several steps whose sequential order must be observed
▪	Action consisting of one or several steps whose sequential order is irrelevant
p. 5	List item
RS-232	Cross-reference to page 5
	Labeling of an operating element on the product (example: socket of the RS-232 interface)
	Warning signs affixed to the product that refer to detailed information in this manual.

1.3 Figures

For better understandability, the colors, proportions, and degree of detail in illustrations can deviate from the actual circumstances. Photographic illustrations may also differ and must not be seen as guaranteed properties.

1.4 Other Applicable Documents

The devices and software tools that are mentioned in this documentation are described in their own manuals.

Product	Document
E-616.SS0/E-616.S0 multi-channel servo controller for piezo tip/tilt mirrors / platforms with SGS	PZ200E User Manual
E-616.SS0G/E-616.S0G multi-channel servo controller for piezo tip/tilt mirrors / platforms with SGS, bench-top	PZ219E User Manual
E-509.S3 sensor/servo controller module	PZ77E User Manual E500T0011 Technical Note for Analog Drivers
E-503.00S piezo amplifier module	PZ62E User Manual
E-501.00 9.5" chassis for modular piezo controller system	
E-518.I3 interface module	E518T0001 Technical Note

1.5 Downloading Manuals

INFORMATION

If a manual is missing or problems occur with downloading:

- Contact our customer service department (p. 35).

INFORMATION

For products that are supplied with software (CD in the scope of delivery), access to the manuals is protected by a password. Protected manuals are only displayed on the website after entering the password.

The password is included on the CD of the product.

For products with CD: Identify the password

1. Insert the product CD into the PC drive.
2. Switch to the Manuals directory on the CD.
3. In the Manuals directory, open the Release News (file including **releasenews** in the file name).
4. Find the user name and the password in the section "User login for software download" in the Release News.

Downloading manuals

1. Open the website **www.pi.ws**.
2. If access to the manuals is protected by a password:
 - a) Click **Login**.
 - b) Log in with the user name and password.
3. Click **Search**.
4. Enter the product code up to the period (e.g., P-882) or the product family (e.g., PICMA® Bender) into the search field.
5. Click **Start search** or press the **Enter** key.
6. Open the corresponding product detail page in the list of search results:
 - a) If necessary: Scroll down the list.
 - b) If necessary: Click **Load more results** at the end of the list.
 - c) Click the corresponding product in the list.
7. Scroll down to the **Downloads** section on the product detail page.

The manuals are displayed under ***Documentation***.

8. Click the desired manual and save it to the hard disk of your PC or to a data storage medium.

2 Safety

In this Chapter

Intended Use.....	5
General Safety Instructions.....	5
Organizational Measures.....	6

2.1 Intended Use

The S-334 is a laboratory device as defined by DIN EN 61010-1. It is intended to be used in interior spaces and in an environment which is free of dirt, oil, and lubricants.

In accordance with its design and realization, the S-334 is intended for the precise positioning and alignment of a mirror in two orthogonal axes with a common pivot point (parallel kinematics). The S-334 can be mounted in any orientation. The S-334 is equipped with strain gauge sensors (SGS).

The intended use of the S-334 is only possible in combination with suitable electronics (p. 14) that is available from PI. The electronics is not included in the scope of delivery of the S-334.

The electronics must provide the required operating voltages. To ensure proper performance of the servo-control system, the electronics must also be able to read out and process the signals from the strain gauge sensors.

2.2 General Safety Instructions

The S-334 is built according to state-of-the-art technology and recognized safety standards. Improper use can result in personal injury and/or damage to the S-334.

- Only use the S-334 for its intended purpose, and only use it if it is in a good working order.
- Read the user manual.
- Immediately eliminate any faults and malfunctions that are likely to affect safety.

The operator is responsible for the correct installation and operation of the S-334.

2.3 Organizational Measures

User manual

- Always keep this user manual available with the S-334.
The latest versions of the user manuals are available for download (p. 3) on our website.
- Add all information from the manufacturer to the user manual, for example supplements or technical notes.
- If you give the S-334 to other users, also include this user manual as well as other relevant information provided by the manufacturer.
- Only use the device on the basis of the complete user manual. Missing information due to an incomplete user manual can result in minor injury and damage to equipment.
- Only install and operate the S-334 after you have read and understood this user manual.

Personnel qualification

The S-334 may only be installed, started up, operated, maintained, and cleaned by authorized and appropriately qualified personnel.

3 Product Description

In this Chapter

Model Overview.....	7
Product View.....	8
Product Labeling.....	10
Scope of Delivery	13
Accessories.....	13
Suitable Electronics.....	14
Control	15
Mirror.....	16

3.1 Model Overview

Two standard versions of the S-334 tip/tilt mirror are available.

Model	Description
S-334.2SD1	High-DynamicsPiezo Tip / Tilt System, 50 mrad, SGS, Sub-D Connector, incl. Mirror
S-334.2SL1	High-DynamicsPiezo Tip / Tilt System, 50 mrad, SGS, LEMO Connectors, incl. Mirror

- For further technical data, see the specifications (p. 37).

3.2 Product View




Figure 1: Schematic product view of the S-334, without transport lock

- 1 Case (material: aluminum)
- 2 Mirror on moving platform
- 3 Cover plate
- 4 Base body (material: titanium)
- 5 Connection cable for piezo actuators and sensors




Figure 2: Schematic diagram of the axes in relation to the cable exit

- 1 Axis 1 (corresponds to tip/tilt axis X on the E-616 controller)
- 2 Axis 2 (corresponds to tip/tilt axis Y on the E-616 controller)
- 3 Cable exit



Figure 3: Maximum displacement in the positive direction of motion around axis 1. At the connected channel 1 of the amplifier, the output voltage U_{piezo} is 100 V. The displacement shown is strongly exaggerated for better understanding.



Figure 4: Maximum displacement in the positive direction of motion around axis 2. At the connected channel 2 of the amplifier, the output voltage U_{piezo} is 100 V. The displacement shown is strongly exaggerated for better understanding.

3.3 Product Labeling




Figure 5: Position of the product labeling, for explanation of the labeling see table









Figure 6: Product labeling on the bottom side of the S-334

The S-334 is labeled as follows:

	Labeling	Description
1	S-334.2Sx1	Product name (example), the characters following the period refer to the model
2	114010244	Serial number (example), individual for each S-334 Meaning of the places (counting from left): 1 = internal information, 2 and 3 = year of manufacture, 4 to 9 = consecutive numbers
3	Country of origin: Germany	Country of origin
4		Warning sign "Observe manual!"
5	WWW.PI.WS	Manufacturer's address (website)
6		Symbol for the protective earth conductor, marks the position of the holes via which the S-334 is to be connected to the protective earth conductor
7		Manufacturer's logo
8		CE conformity mark
9		Old equipment disposal (p. 43)

S-334.2SD1: Labeling of the Sub-D 25 (m) connector




Figure 7: Sub-D 25 (m) connector on the connection cable of the S-334.2SD1

1



Warning sign "Residual voltage": Notice of risk of electric shock (p. 5)

S-334.2SL1: Labeling of the connection cables




Figure 8: The cables of the S-334.2SL1

- 1 Sensor connection of axis 1, labeled **Sensor 1**
- 2 Sensor connection of axis 2, labeled **Sensor 2**
- 3 Piezo connection of axis 1, labeled **PZT 1**
- 4 Piezo connection of axis 2, labeled **PZT 2**
- 5 Piezo connection for 100 V fixed voltage, labeled **PZT 3**

3.4 Scope of Delivery

Item ID	Components
S-334	Piezo tip/tilt mirror according to order
PZ255E	User manual (this document) in printed form
S334T0001	Technical Note with handling information in printed form
S334E0002	Transport lock

3.5 Accessories

Order number	Description
Only S-334.xxL:	
E-518.13	Interface Module, 3 Channels, TCP/IP, USB and RS-232 Interfaces
P-891.01	Extension Cable for Piezo Voltage, LEMO Connectors, 1 m
P-891.02	Extension Cable for Piezo Voltage, LEMO Connectors, 2 m
P-891.03	Extension Cable for Piezo Voltage, LEMO Connectors, 3 m
P-891.05	Extension Cable for Piezo Voltage, LEMO Connectors, 5 m
P-891.10	Extension Cable for Piezo Voltage, LEMO Connectors, 10 m
P-892.01	Extension Cable, for Strain Gauge Sensors, LEMO Connectors, 1 m
P-892.02	Extension Cable, for Strain Gauge Sensors, LEMO Connectors, 2 m
P-892.03	Extension Cable, for Strain Gauge Sensors, LEMO Connectors, 3 m
P-892.05	Extension Cable, for Strain Gauge Sensors, LEMO Connectors, 5 m
P-892.10	Extension Cable, for Strain Gauge Sensors, LEMO Connectors, 10 m

➤ To order, contact our customer service department (p. 35).

3.6 Suitable Electronics

Tip/Tilt Mirror	Controller	Amplifier	Case	Display
S-334.2SD1	E-616.SS0 Multi-Channel Servo Controller for Piezo Tip/Tilt Mirror Platforms with SGS and Differential Drive		-	-
	E-616.SS0G Multi-Channel Servo Controller for Piezo Tip/Tilt Mirror Platforms with SGS and Differential Drive, Benchtop Device		205 mm × 105 mm × 54.1 mm	-
S-334.2SL1	E-509.S3 Sensor / Piezo Servo-Control Module, SGS Sensors, 3 Channels	E-503.00S Piezo Amplifier Module, -30 to 130 V, 2 Channels, Modified E-503.00 for S-330, S-334, S-340 Tip/Tilt Systems, with One Fixed Voltage of +100V, Two Variable Voltages	E-501.00 9½" Chassis for Modular Piezo Controller System, 1 to 3 Channels	Optional: E-518.I3 Interface Module, 3 Channels, TCP/IP, USB and RS-232 Interfaces

➤ To order, contact our customer service department (p. 35).

3.7 Control




Figure 9: Differential drive of the piezo tip/tilt system, functional principle using the tilting of a single axis as an example

- 1 Piezo actuator 1 of the axis
- 2 Piezo actuator 2 of the axis
- 3 Moving platform with mirror

The S-334 is a piezo tip/tilt mirror with a differential piezo drive. Tip/tilt motions in two axes are realized by the pairwise interconnection of a total of four piezo actuators.

Both pairs of actuators are electrically connected in such a way that when the piezo voltage U_{Piezo} is changed, an increased voltage is applied to one actuator of a pair while the applied voltage on the other actuator is decreased by the same amount. The actuator with the increased voltage expands while the other actuator with the decreased voltage contracts. This produces the tip/tilt motion.

For a simplified representation of the functional principle, only one axis is shown in the figure above. The moving platform is shown with a tilt of 0° .

When the control input voltage U_{In} increases, piezo actuator 1 expands and piezo actuator 2 contracts. This produces a tilt in the positive direction.

Due to the type of interconnection, both actuators of a pair always move in opposite directions. It is thus impossible to command linear motions in the Z axis.

The position of the Z axis can change with temperature fluctuations, however: Due to the symmetrical design of the piezo tip/tilt mirror, temperature fluctuations do not cause the moving platform to tilt but cause the length of the piezo actuators to change evenly in the direction of the Z axis.

Most applications are not very sensitive to such deviations, as long as the tip/tilt angle does not change.

The four piezo actuators of the S-334 are each equipped with a strain gauge sensor. For this reason, a servo loop with a sensor channel must be available to each pair of actuators, in addition to the amplifier channel.

3.8 Mirror

The S-334 is equipped with a factory-mounted mirror with the following characteristics:

- Diameter 10 mm
- 2 mm thickness
- Surface accuracy $\lambda/10$
- Surface quality 20-10
- Parallelism 30 arc seconds

4 Unpacking

INFORMATION

The S-334 is delivered with a transport lock.

- To avoid scratches on the mirror surface during mounting, it is recommended not to remove the transport lock until after mounting (p. 23).

1. Unpack the S-334 with care.
2. Compare the contents against the items covered by the contract and against the packing list.
3. Inspect the contents for signs of damage. If there is any sign of damage or missing parts, contact PI immediately.
4. Keep all packaging materials in case the product needs to be returned.

5 Installation

In this Chapter

General Notes on Installation	19
Mounting the S-334 and Connecting it to a Protective Earth Conductor.....	21
Removing the Transport Lock	23
Connecting the S-334 to the Controller	24

5.1 General Notes on Installation

CAUTION



Dangerous voltage and residual charge on piezo actuators!

The S-334 is driven by piezo actuators. Temperature changes and compressive stresses can induce charges in piezo actuators. After being disconnected from the electronics, piezo actuators can also stay charged for several hours. Touching or short-circuiting the contacts in the connector of the S-334 can lead to minor injuries. In addition, the piezo actuators can be destroyed by an abrupt contraction.

- Do **not** open the S-334.
- Discharge the piezo actuators of the tip/tilt mirror before installation:
Connect the tip/tilt mirror to the switched-off PI controller, which is equipped with an internal discharge resistor.
- Do **not** pull the connector out of the electronics during operation.



In the case of tip/tilt mirrors with Sub-D connectors:

Touching the contacts in the connector can lead to an electric shock (max. 100 V DC) and minor injuries.

- Do **not** touch the contacts in the connector.
- Secure the connector of the tip/tilt mirror with screws against being pulled out of the controller.

NOTICE**Destruction of the piezo actuator due to electric flashovers!**

Using the S-334 in environments that increase the electrical conductivity can lead to the destruction of the piezo actuator by electric flashovers. Electric flashovers can be caused by moisture, high humidity, liquids, and conductive materials (e.g., metal dust). In addition, electric flashovers can also occur in certain air pressure ranges due to the increased conductivity of the air.

- Avoid operating the S-334 in environments that can increase the electric conductivity.
- Only operate the S-334 within the permissible ambient conditions and classifications (p. 39).

NOTICE**Destruction of the piezo actuator due to short-circuiting without a discharge resistor!**

When a charged piezo actuator is short-circuited without a discharge resistor, this can lead to a contraction shock and thus to the destruction of the piezo ceramic.

- Only discharge the S-334 according to the instructions in "Discharging S-334" (p. 28).

NOTICE**Damage when the mirror is removed!**

The mirror of the S-334 may only be replaced by PI. Otherwise, the S-334 can be damaged.

- Do **not** remove the mirror of the S-334.
- If you need a different mirror, contact our customer service department (p. 35).

NOTICE**Warping of the S-334 due to mounting on uneven surfaces!**

Mounting the S-334 on an uneven surface can warp the S-334. Warping reduces the accuracy.

- Mount the S-334 on an even surface. The recommended evenness of the surface is $\leq 30 \mu\text{m}$.
- For applications with large temperature changes:
Only mount the S-334 on surfaces that have the same or similar thermal expansion properties as the S-334.

NOTICE**Damage due to unsuitable cables!**

Unsuitable cables can damage the S-334 and the electronics.

- Only use cables provided by PI for connecting the S-334 to the electronics.

5.2 Mounting the S-334 and Connecting it to a Protective Earth Conductor

The S-334 can be mounted as follows:

- Mounting on a surface with two screws
- Clamping in a clamping holder



Figure 10: M3 holes of the S-334.2Sx1 for mounting on a surface and/or protective earth connection



Figure 11: Clamping width for mounting the S-334.2Sx1 by clamping

NOTICE



Damage due to incorrect mounting!

Incorrect mounting can damage the S-334.

- When mounting the S-334 by clamping, observe the maximum clamping width of 20 mm and the maximum clamping force of 30 N.
- Only mount the S-334 according to the instructions in this manual.

INFORMATION

The contact of the S-334 with the protective earth conductor is made as follows:

- Two M3 holes in the base body of the S-334, marked with the symbol for the protective earth conductor (⏚)
- Suitable conductive M3 screws

When the S-334 is mounted on a surface via the two M3 holes:

- The surface is sufficiently conductive and connected to the protective earth conductor.

When the S-334 is mounted by clamping:

- The protective earth conductor is connected to both M3 holes of the S-334 via two screws.

INFORMATION

- Observe the applicable standards for connecting the protective earth conductor.

Requirements

- ✓ You have read and understood the general notes on installation (p. 19).
- ✓ The S-334 is **not** connected to the controller.
- ✓ You have accounted for the space required for a cable routing free of kinks and in accordance with regulations.

Tools and accessories

- For the dimensions of the tip/tilt mirror and the position and depth of the M3 holes, see "Dimensions" (p. 40).
- When the S-334 is mounted on a surface: You have provided a suitable surface:
 - The surface must be connected to a protective earth conductor.
 - Two through holes for M3 screws are present.
 - The holes for accommodating the screws have to be sufficiently conductive to ensure the proper functioning of the protective earth conductor.
 - The contact resistance at all connection points relevant for mounting the protective earth conductor is $<0.1 \Omega$ at 25 A.
 - The evenness of the surface is $\leq 30 \mu\text{m}$.
- When the S-334 is mounted by clamping:
 - Suitable clamping holder for optical components, permissible clamping width of the S-334 see figure above
- Suitable protective earth conductor: Cross-sectional area of the cable $\geq 0.75 \text{ mm}^2$
- Two electrically conductive M3 screws of appropriate length (p. 40)
- Suitable tools

Mounting the S-334 on a surface and connecting it to a protective earth conductor

1. Align the S-334 on the surface so that the M3 holes in the S-334 and the surface overlap.
2. Introduce the two screws through the holes in the surface into the base body of the S-334 from below.
3. Tighten the two screws.
 - Maximum screw-in depth: 4 mm
 - Maximum torque: 1.1 Nm

4. Check that the S-334 fits without backlash.

Mounting the S-334 with a clamping holder and connecting it to the protective earth conductor

- Mount the S-334 in the clamping holder.
 - Permissible clamping width, see figure above
 - Maximum clamping force: 30 N
- Connect the protective earth conductor to the M3 holes marked with the symbol for the protective earth conductor (⏚) in the base body of the S-334 via two screws.
 - Maximum screw-in depth: 4 mm
 - Maximum torque: 1.1 Nm

5.3 Removing the Transport Lock

The S-334 is delivered with a transport lock.



Figure 12: Transport lock of the S-334, the arrows mark the recommended adhesion points

Tools and accessories

- Suitable tool for removing the transport lock
- For removing the adhesion points, if necessary: isopropanol and cotton swab

Requirements

- ✓ You have mounted the S-334 (p. 21).

Removing the transport lock

1. Carefully remove the transport lock.
2. If necessary: Remove the adhesive strips for the transport lock and the adhesive residue with a cotton swab soaked in isopropanol.
3. Keep the transport lock in case the product needs to be transported again later.

5.4 Connecting the S-334 to the Controller

INFORMATION

Systems consisting of an S-334 and controller are calibrated at the factory to achieve optimum performance.

- Note the assignment of the axes to the controller channels, which is specified on the calibration label of the piezo servo controller.

Requirements

- ✓ You have read and understood the general notes on installation (p. 19).
- ✓ You have installed a suitable controller (p. 14).
- ✓ You have read and understood the user manual of the controller.
- ✓ The controller is switched off.

Connecting the S-334.2SD1 to the E-616 controller

1. Plug the connector of the S-334.2SD1 into the corresponding socket of the controller (see user manual of the controller).
2. Use the integrated screws to secure the connection against accidental disconnection.

Connecting the S-334.2SL1 to E-50x modules

1. Connect the piezo connections of the S-334.2SL1 to the E-503.00S piezo amplifier module as follows:
 - **PZT 1** to **PZT** for channel 1 (**CH1**)
 - **PZT 2** to **PZT** for channel 2 (**CH2**)
 - **PZT 3** to **PZT** for channel 3 (**CH3**)
2. Connect the sensor connections of the S-334.2SL1 to the E-509.S3 servo controller module as follows:
 - **Sensor 1** to **SENSOR** for channel 1 (**SERVO 1**)
 - **Sensor 2** to **SENSOR** for channel 2 (**SERVO 2**)

6 Start-Up and Operation

In this Chapter

General Notes on Start-Up and Operation	25
Operating the S-334	27
Discharging the S-334	28

6.1 General Notes on Start-Up and Operation

CAUTION



Risk of electric shock if the protective earth conductor is not connected!

If a protective earth conductor is not or not properly connected, dangerous touch voltages can occur on the S-334 in the case of malfunction or failure of the system. If touch voltages exist, touching the S-334 can result in minor injuries from electric shock.

- Before start-up, establish contact between the S-334 and the protective earth conductor (p. 21).
- Do **not** remove the protective earth conductor during operation.
- Make sure that the contact resistance at all connection points relevant for mounting the protective earth conductor is $<0.1 \Omega$ at 25 A.
- If the protective earth conductor has to be temporarily removed (e.g., for modifications), reconnect the S-334 to the protective earth conductor before starting it up again.

NOTICE



Destruction of the piezo actuator due to electric flashovers!

Using the S-334 in environments that increase the electrical conductivity can lead to the destruction of the piezo actuator by electric flashovers. Electric flashovers can be caused by moisture, high humidity, liquids, and conductive materials (e.g., metal dust). In addition, electric flashovers can also occur in certain air pressure ranges due to the increased conductivity of the air.

- Avoid operating the S-334 in environments that can increase the electric conductivity.
- Only operate the S-334 within the permissible ambient conditions and classifications (p. 39).

NOTICE**Reduced lifetime of the piezo actuator due to permanently high voltage!**

The permanent application of a high static voltage to piezo actuators leads to a considerable reduction in the lifetime of the piezo ceramic of the S-334.

- When the S-334 is not used but the controller remains switched on to ensure temperature stability, discharge the S-334 (p. 28).

NOTICE**Impermissible mechanical load due to unsuitable control!**

The S-334 is designed for the following control:

- Maximum range of piezo voltage: 0 to 100 V
- Minimum rise time of the piezo voltage: 10 ms
- Maximum frequency with sinusoidal control signal: 50 Hz
- Maximum frequency with 10 ms rise time of the piezo voltage: 30 Hz

An unsuitable control of the S-334 results in an impermissibly high mechanical load and can cause damage to the S-334.

- Set the control signal so that the piezo voltage stays within the limit values.
- Only use the S-334 in closed-loop operation.

NOTICE**Damage from start-up with transport lock!**

If the transport lock has **not** been removed before start-up, the mirror can collide with the transport lock. Collisions can cause damage to the mirror.

- Remove the transport lock before you start up the S-334 (p. 23).

INFORMATION

Systems consisting of an S-334 and controller are calibrated at the factory to achieve optimum performance.

- Note the assignment of the axes to the controller channels, which is given by the calibration label of the piezo servo controller.
- Do **not** change the dynamic servo-control parameters of the controller (e.g. P term and I term) for a system calibrated by PI.

Only when replacing system components:

- If the controller or the S-334 has to be replaced, recalibrate the axis displacement (see controller manual) or contact our customer service department (p. 35).
- If the controller has to be replaced, set the dynamic servo-control parameters (see controller manual).

INFORMATION

Sound and vibration (e.g., footfall, impacts) can be transmitted to the S-334 and can affect its performance with regard to position stability.

- Avoid sound and vibration during operation of the S-334.

INFORMATION

In the case of a S-334 tip/tilt mirror with Sub-D connector, ground loops can occur when the tip/tilt mirror is grounded via its protective earth connection as well as by the shield of the connection cable for the electronics.

- If a ground loop occurs, contact our customer service department (p. 35).

INFORMATION

Optionally, the S-334.2SL1 can also be operated from a PC via a program such as PIMikroMove using the digital E-518.I3 interface module.

INFORMATION

The expansion of the piezo actuators depends on the ambient temperature and can vary by up to 10 % in the given temperature ranges (p. 39).

6.2 Operating the S-334

Prerequisites

- ✓ You have read and understood the general notes on start-up and operation (p. 25).
- ✓ You have read and understood the user manual of the controller.
- ✓ If you use the E-518.I3 interface module: You have read and understood the user manual of the interface module and the PC software.
- ✓ You have properly installed the tip/tilt mirror (p. 19).
- ✓ The controller and the required PC software have been installed. All connections with the controller have been set up (see user manual of the controller).
- ✓ You have removed the transport lock (p. 23).
- ✓ You have ensured that the servo mode is switched on for all axes of the S-334 on the controller.

Operating the S-334

- For start-up and operation of the S-334, follow the instructions in the manual of the controller used.

6.3 Discharging the S-334

The S-334 must be discharged in the following cases:

- Before installation
- If the S-334 is not used but the controller remains switched on to ensure temperature stability
- Before pulling out the connector (e.g., before cleaning and transport of the S-334 and for modifications in the application)

The S-334 is discharged through the internal discharge resistor of the controller from PI.

Discharging a S-334 that is connected to the controller

1. Switch the controller off.
2. Wait at least 10 seconds before pulling out the connector.

Discharging a S-334 that is not connected to the controller

- Connect the S-334 to the switched-off controller from PI.

7 Maintenance

In this Chapter

General Notes on Maintenance.....	29
Cleaning the S-334	30
Packing the S-334 for Transport	30

7.1 General Notes on Maintenance

NOTICE



Misalignment due to loosening screws!

The S-334 is maintenance-free and achieves its positioning accuracy as a result of the optimum alignment of mechanical components and piezo actuators. Loosened screws cause a loss in positioning accuracy.

- Only loosen screws according to the instructions in this manual.
- Do **not** open the S-334.

NOTICE



Damage when the mirror is removed!

The mirror of the S-334 may only be replaced by PI. Otherwise, the S-334 can be damaged.

- Do **not** remove the mirror of the S-334.
- If you need a different mirror, contact our customer service department (p. 35).

7.2 Cleaning the S-334

NOTICE

**Damage due to incorrect cleaning!**

The mirror of the S-334 can be damaged from applying force during cleaning.

- Only clean the mirror when actually necessary.
- Avoid exerting any force on the mirror during cleaning.
- Do **not** use compressed air.

Requirements

- ✓ You have discharged the piezo actuators of the S-334 (p. 28).
- ✓ You have disconnected the S-334 from the controller.

Tools and accessories

- Bellows
- Optic brush

Cleaning the S-334

- Clean the surfaces and the mirror of the S-334 with bellows and/or an optic brush without exerting force.
- Do **not** use compressed air for cleaning.
- Do **not** do any ultrasonic cleaning.

7.3 Packing the S-334 for Transport

Accessories

- Easy to remove, double-sided adhesive tape
- Suitable foil
- Transport lock (p. 17)

Packing the S-334 for transport

Figure 13: Transport lock of the S-334, the arrows mark the recommended adhesion points

1. Fix the transport lock on the S-334 with double-sided adhesive tape as shown in the figure above.



Figure 14: S-334 in the original packaging, with dessiccants

2. Pack the S-334 as shown in the figure above, if possible.

8 Troubleshooting

Problem	Possible Causes	Solution
No or uncontrolled motion	<ul style="list-style-type: none"> ▪ Cable not connected correctly ▪ Controller defective ▪ Cable defective ▪ Piezo ceramic defective after electric flashover 	<ul style="list-style-type: none"> ➤ Check the cable connections (p. 24). ➤ Contact our customer service department (p. 35).
	<ul style="list-style-type: none"> ▪ S-334 damaged due to open-loop operation 	<ul style="list-style-type: none"> ➤ Contact our customer service department (p. 35).
Reduced accuracy	Warped base body	Only mount the S-334 on surfaces with the following characteristics: <ul style="list-style-type: none"> – Evenness of at least 30 µm – The thermal expansion properties are similar to those of the S-334 (e.g., surfaces made of titanium)
	Impermissible clamping of the S-334	<ul style="list-style-type: none"> ➤ Observe the permissible clamping width and clamping force (p. 21).
	S-334 or controller has been replaced	<p>After the S-334 or controller has been replaced, it is necessary to recalibrate the axis displacement.</p> <ul style="list-style-type: none"> ➤ Perform a recalibration of the axis displacement (see controller manual) or contact our customer service department (p. 35).
	Only with S-334.2SL1: Axes were mixed up during connection	<ul style="list-style-type: none"> ➤ Note the assignment of the axes to the controller channels, which is given by the calibration label of the piezo servo controller.
	Operating temperature outside of the permissible range (p. 37)	<ul style="list-style-type: none"> ➤ Contact our customer service department (p. 35).

If the problem that occurred with your system is not listed in the table above or cannot be solved as described, contact our customer service department (p. 35).

9 Customer Service

For inquiries and orders, contact your PI sales engineer or send us an email (info@pi.ws).

- If you have questions concerning your system, have the following information ready:
 - Product codes and serial numbers of all products in the system
 - Firmware version of the controller (if present)
 - Version of the driver or the software (if present)
 - Operating system on the PC (if present)
- If possible: Take photographs or make videos of your system that can be sent to our customer service department if requested.

The latest versions of the user manuals are available for download (p. 3) on our website.

10 Technical Data

In this Chapter

Specifications	37
Dimensions	40
Pin Assignment	41

10.1 Specifications




10.1.1 Data Table

	S-334.2SL1 / S-334.2SD1	Unit	Tolerance
Active axes	θ_x, θ_y		
Motion and positioning			
Integrated sensor	SGS		
Closed-loop tilt angle	50	mrاد	
Open-loop resolution	0.5	μrad	typ.
Closed-loop resolution	5	μrad	typ.
Linearity	0.05	%	typ.
Repeatability	5	μrad	typ.
Mechanical properties			
Resonant frequency under load (with standard mirrors)	0.7	kHz	±20 %
Load capacity	0.2	N	max.
Distance of pivot point to platform surface	4	mm	±0.1 mm
Drive properties			
Ceramic type	PICMA® P-885		
Electrical capacitance per axis	3	μF	±20 %
Miscellaneous			
Operating temperature range	0 to 50	°C	

	S-334.2SL1 / S-334.2SD1	Unit	Tolerance
Material case	Titanium		
Mass	0.065	kg	±5 %
Cable length	2	m	±0.1 m
Sensor / voltage connection	LEMO connectors / Sub-D connector 25-pin		
Recommended controller / amplifier	Version with Sub-D connector: E-616 servo controller for tip/tilt mirror systems Version with LEMO connectors: E-500 modular piezo controller system with E-503.00S amplifier module (three channels) and E-509 controller		

10.1.2 Maximum Ratings

S-334 tip/tilt mirrors are designed for the following operating data:

Maximum Operating Voltage	Maximum Operating Frequency ¹ (with Standard Mirror)	Maximum Power Consumption ²
		
0 to +100 V	30 Hz ³	11 W / axis

¹ To ensure stable operation, the maximum operating frequency has been defined as around one third of the mechanical resonant frequency.

² The heat that is generated by the piezo actuator during dynamic operation limits the value for maximum power consumption.

Details can be found at the following website:

<http://piceramic.com/piezo-technology/properties-piezo-actuators/electrical-operation.html>

³ With a rise time of the piezo voltage of 10 ms. Maximum operating frequency (with standard mirror) for sinusoidal control signal: 50 Hz

10.1.3 Ambient Conditions and Classifications

The following ambient conditions and classifications must be observed for the S-334:

Area of application	For indoor use only
Maximum altitude	2000 m
Air pressure	1100 hPa to 0.1 hPa (equivalent of roughly 825 Torr to 0.075 Torr)
Relative humidity	Highest relative humidity 80 % for temperatures up to 31 °C Decreasing linearly to 50 % relative humidity at 40 °C
Storage temperature	–20 °C to 80 °C
Transport temperature	–25 °C to 85 °C
Overvoltage category	II
Protection class	I
Degree of pollution	1
Degree of protection according to IEC 60529	IP20

10.2 Dimensions

Dimensions in mm. Note that the decimal places are separated by a comma in the drawings.




Figure 15: The dimensions for S-334.2SD1 and S-334.2SL1 are identical

L cable length




Figure 16: S-334.2SL1: Position and dimensions of the cable splitter box

10.3 Pin Assignment

10.3.1 S-334.2SD1: Sub-D 25 (m) Piezo and Sensor Connection




Figure 17: Piezo and sensor connector, Sub-D 25 (m), contact side

Pin	Signal	Function
1	-	-
2	-	-
3	SGS Y+	SGS signal axis 2 (positive)
4	SGS Ref Y	SGS reference axis 2
5	SGS X+	SGS signal axis 1 (positive)
6	SGS Ref X	SGS reference axis 1
7	-	-
8	-	-
9	-	-
10	GND	Ground
11	PZT X	Piezo voltage axis 1
12	PZT Y	Piezo voltage axis 2
13	100 V fix	100 V fixed voltage
14	-	-
15	-	-
16	SGS Y-	SGS signal axis 2 (negative)
17	SGS GND Y	Ground SGS signal axis 2
18	SGS X-	SGS signal axis 1 (negative)
19	SGS GND X	Ground SGS signal axis 1
20	-	-

Pin	Signal	Function
21	-	-
22	-	-
23	PZT GND	Ground piezo voltage
24	PZT GND	Ground piezo voltage
25	PZT GND	Ground piezo voltage

10.3.2 S-334.2SL1: LEMO Piezo and Sensor Connections



Figure 18: Sensor connection: LEMO connector FFA.0S.304.CLAC32Y, contact side

Pin	Signal	Function
1	SGS Ref	SGS reference
2	SGS-	SGS signal (negative)
3	SGS+	SGS signal (positive)
4	SGS GND	Ground SGS signal

PZT




Figure 19: Piezo connector

Signal	Function	Connector Shell
PZT	Piezo Voltage	Ground

11 Old Equipment Disposal

In accordance with EU law, electrical and electronic equipment may not be disposed of in EU member states via the municipal residual waste.

Dispose of your old equipment according to international, national, and local rules and regulations.

In order to fulfil its responsibility as the product manufacturer, Physik Instrumente (PI) GmbH & Co. KG undertakes environmentally correct disposal of all old PI equipment made available on the market after 13 August 2005 without charge.

Any old PI equipment can be sent free of charge to the following address:

Physik Instrumente (PI) GmbH & Co. KG
Auf der Roemerstr. 1
D-76228 Karlsruhe, Germany



12 EC Declaration of Conformity

For the S-334, an EC Declaration of Conformity has been issued in accordance with the following European directives:

Low Voltage Directive

EMC Directive

RoHS Directive

The applied standards certifying the conformity are listed below.

Safety (Low Voltage Directive): EN 61010-1

EMC: EN 61326-1

RoHS: EN 50581

