
Q-Motion® Stages

MP139EK Short Instructions
Q-5xx / Q-6xx

Version: 3.0.0 Date: 23.03.2017

User Information

These short instructions contain an overview of the most important safety instructions for stages with Q-Motion® drives with the product codes given above (x: arbitrary number).

Downloading and Reading the Manual

The actions during installation, start-up and maintenance require additional information from the manuals of the stage and/or the controller used.

Manuals may be titled as follows: „User Manual“, „Technical Note“.

Downloading manuals from the website:

1. Make a note of the product code on the product (e.g. Q-522.140).
2. Open the website www.pi.ws.
3. For password-protected manuals (CD in the scope of delivery of the product):
 - a. Insert the CD of the product into the PC.
 - b. Make a note of the user name and password in the **Releasenews** file in the **\Manuals** folder on the CD.
 - c. Click **Login**.
 - d. Log in with the user name and password.
4. Click **Search**.
5. Enter the product code up to the period (e. g., Q-522) into the search field.
6. Click **Start search** or press the **Enter** key.
7. If necessary: Click **Load more results** at the end of the page.
8. Click on the corresponding product to open the product details page.
9. Scroll down to the **Download** section.
10. Click the desired manual in the **Documentation** subsection and save it onto the hard disk of the PC or onto a data storage medium.

If you cannot find the manual you are looking for, or if you have any questions: Contact our customer service department via info@pi.ws.

PI miCos GmbH, Freiburger Strasse 30, 79427 Eschbach, Germany
Phone +49 7634 5057-0, Fax +49 7634 5057-99, Email info@pimicos.com, www.pi.ws

MP139EK 3.0.0 03/2017.0.5

Safety Instructions

Intended Use

The stage is a laboratory device as defined by DIN EN 61010. It is intended to be used in interior spaces and in an environment which is free of dirt, oil and lubricants.

In accordance with its design, the stage is intended for positioning and adjusting loads at different velocities in interval operation. The stage is not intended for applications in areas in which a failure would present severe risks to human beings or the environment.

The intended use of the stage is only possible when completely mounted and connected and only in combination with a suitable controller. The stage may only be installed, operated, maintained and cleaned by authorized and appropriately qualified personnel.

Installation

Electrostatic discharge can damage the stage.

- ▶ If the stage is supplied with an ESD protective cap: Only remove the ESD protective cap from the connector when connecting the stage to the controller.

Mechanical forces can damage or misalign the stage.

- ▶ Avoid shocks and drops.
- ▶ Observe the maximum permissible forces (see manual).
 - ▶ **Linear stage:** When the motion axis is oriented vertically, the load must be lower than the holding force of the drive (see manual).
 - ▶ **Rotary stage:** When the stage is mounted vertically, the load must be lower than the maximum torque of the drive (see manual).
- ▶ Include the masses of the moved stages in multi-axis systems in the calculations.

Dirt, oil, lubricants, and condensation will render both the drive and the position sensor of the stages inoperable.

- ▶ Keep the piezomotor free from lubricants.
- ▶ Do **not** touch the scale of the encoder.
- ▶ Keep the stage free from dirt and condensation.

If the scale of the encoder is accessible, it can be damaged by sharp objects.

- ▶ Treat the scale of the encoder with extreme caution.

Incorrectly fitted screws and locating pins can cause damage.

- ▶ Observe the maximum torque for mounting (see manual).
- ▶ Do **not** allow screw heads to protrude.
- ▶ Select the screw length according to the depth of the mounting holes.
- ▶ Do not insert locating pins too deeply into the motion platform.

Collisions can damage the stage, the load to be moved and the environment.

- ▶ Mount the stage and the load so that the load cannot get jammed or blocked, or collide with objects in the workspace.

Unsuitable mounting can warp the stage and reduce the accuracy.

- ▶ Mount the stage on an even surface with similar thermal expansion properties (for recommended evenness, see manual).

Unsuitable cables can cause damage to the electronics and can affect the performance of the stage.

- ▶ Only use original parts from PI (see manual).

Heat produced during operation can affect your application.

- ▶ Install the stage so that the application is not affected by the dissipated heat.
- ▶ Ensure sufficient ventilation and heat dissipation from contact surfaces.

Start-Up

If a protective earth conductor is not or not properly connected, touching the stage can result in minor injuries from electric shock in the case of a malfunction.

- ▶ Only operate the stage with a properly connected protective earth conductor.
- ▶ Do **not** remove the protective earth conductor from the stage during operation.
- ▶ Observe the applicable standards for mounting the protective earth conductor.

MP139EK 3.0.0 03/2017.0.5

The surface of the stage and the surrounding area can heat up during operation. Touching the stage can lead to minor injuries from burning. Overheating can damage the stage.

- ▶ Cool the stage and the surrounding area or install touch protection.
- ▶ Select the operating time according to the operating frequency in step mode (see manual).

Moisture, liquids and electrically conductive materials (e.g., metal dust) that penetrate the stage can destroy the piezomotor of the stage.

- ▶ Only operate the stage within the permissible ambient conditions (see manual).
- ▶ Do **not** operate a vacuum-compatible stage during evacuation.

Operating voltages that are too high or incorrectly connected can cause damage to the stage.

- ▶ Only use compatible electronics.
- ▶ Observe the operating voltage range of the stage (see manual).
- ▶ Observe the correct pin assignment (see manual).

Continuously high static voltages can reduce the lifetime of the piezomotor.

- ▶ Switch off the electronics if the stage is not in use.

A high operating frequency in step mode and continuous use of a small working range can cause considerable wear of the stage.

- ▶ Reduce the operating frequency and duration.
- ▶ Stop the motion immediately if the electronics malfunction.
- ▶ If possible: Select another part of the travel range as working range at regular intervals.

Handling of Vacuum-Compatible Products

If the stage is vacuum-compatible, attention must be paid to cleanliness.

- ▶ Only touch the stage with powder-free gloves.
- ▶ If necessary, wipe the stage clean.

Mounting the Stage and Connecting it to the Protective Earth Conductor


NOTICE!

- ▶ Damage from incorrect mounting.
 - ▶ Observe the safety instructions in the „Installation“ section.

NOTICE!

- ▶ Damage from mechanical actions or dirt.
 - ▶ Do not touch the ceramic rail of the Q-Motion® drive.
 - ▶ Treat the scale of the encoder with extreme caution.

Scale and ceramic rail of the Q-Motion® drive

MP139EK 3.0.0 03/2017.0.5

Suitable protective earth conductor

- ▶ Cross-section of the protective earth conductor $\geq 0.75 \text{ mm}^2$.

Mounting the stage

- ▶ To fix the stage to the surface, use sufficiently conductive screws in all mounting holes.
- ▶ Tighten the mounting screws with at least three rotations (torque see manual).
- ▶ Secure the screws against unintentional loosening, e. g., with locking paint.


Checking the connection of the protective earth conductor

- ▶ Make sure that the contact resistance is $< 0.1 \Omega$ at 25 A at all protective earth conductor connections.
- ▶ If ground loops occur, contact the PI customer service department.


Affixing the Load




- ▶ Affix the load to the mounting holes provided.

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