Industrial Applications

LMP - Laser Micromachining

Many laser micromachining processes used, for example, in the production of photovoltaic cells, electronic devices, and medical microelements require a small focal point to achieve high intensity and produce fine features. The V-308 enables high quality laser micromachining by providing nanometer precision control of the focal height. Even with irregular workpiece surfaces the V-308 allows fast compensation of the focusing optics according to the required focus height to avoid loss of quality. The compensation can be achieved by precharacterization or by using a high-speed autofocus sensor. This allows high volume production in processes where a narrow depth of focus is required.





LMP – Laser Ablation

The laser ablation of the surface of engineering ceramics, such as silicon or sapphire, requires a narrow depth of focus to ensure a small focus point that leads to an accurate control of material removal and limits damage. The laser spot defines the fluence, i.e., the energy delivered per unit area, (the energy delivered per unit area) and its regulation is needed to control the ablation threshold. When material is removed, the distance between the surface being processed and the current focal point change. With the V-308, it is possible to control exactly how deep in the material the process takes place by adjusting the focal point. In addition, the processing efficiency and the quality of the sidewalls within the features can be optimized.

Wafer Inspection

As semiconductor devices move towards a higher integration and density with the advances in semiconductor materials and processes, the critical defect size decreases. There are two factors that affect the topology of the wafer surface: The surface tension induced by IC manufacturing and the surface roughness. To guickly obtain a high-quality surface image, defect inspection tools must operate at high resolutions and high inspection speeds. A precise and fast and fast autofocus ability is critical. The V-308 allows a fast and precise compensation to maintain depth of focus. It has the added advantage of an extended travel range, eliminating the need for an additional coarse moving stage.



Headquarters

GERMANY

Physik Instrumente (PI) , GmbH & Co. KG Auf der Roemerstrasse 1 76228 Karlsruhe Phone +49 721 4846-0 Fax +49 721 4846-1019 info@pi.ws www.pi.ws

Pl miCos GmbH Freiburger Strasse 30 79427 Eschbach Phone +49 7634 5057-0 Fax +49 7634 5057-99 info@pimicos.com www.pi.ws

PI Ceramic GmbH

Lindenstrasse 07589 Lederhose Phone +49 36604 882-0 Fax +49 36604 882-4109 info@piceramic.com www.piceramic.com



© Physik Instrumente (PI) GmbH & Co. KG

All contents, including texts, graphics, data etc., as well as their layout, are subject to copyright and other protective laws. Any copying, modification or redistribution in whole or in parts is subject to a written permission of Pl.

ISRAEL

POB 984

1 Hataasia St.

ACS Motion Control Ltd.

Migdal HaEmek, 2307037

Phone +972-4-6546440

Fax +972-4-6546443

Ramat Gabriel Industrial Park

info@acsmotioncontrol.com

www.acsmotioncontrol.com

Although the information in this document has been compiled with the greatest care, errors cannot be ruled out completely. Therefore, we cannot guarantee for the information being complete, correct and up to date. Illustrations may differ from the original and are not binding. PI reserves the right to supplement or change the information provided without prior notice.



WWW.PI.WS

WWW.PI.WS





ACS Motion Control PI Subsidiaries

USA (East) & CANADA	USA (West) & MEXICO
PI (Physik Instrumente) L.P.	PI (Physik Instrumente) L.P.
Auburn, MA 01501	Irvine, CA 92620
www.pi-usa.us	www.pi-usa.us
USA (San Francisco Bay Area)	UK & IRELAND
PI (Physik Instrumente) L.P.	PI (Physik Instrumente) Ltd.
Sausalito, CA 94965	Cranfield, Bedford
www.pi-usa.us	www.physikinstrumente.co.uk
ITALY	NETHERLANDS
Physik Instrumente (PI) S. r. l.	PI Benelux B.V.
Bresso	Sint-Oedenrode
www.pionline.it	www.pi.ws/benelux
FRANCE	SPAIN
PI France SAS	Micos Iberia S.L.
Aix-en-Provence	Vilanova i la Geltrú
www.pi.ws	www.pimicos.es
JAPAN	
PI Japan Co., Ltd.	PI Japan Co., Ltd.
Tokyo	Osaka
www.pi-japan.jp	www.pi-japan.jp
CHINA	
Physik Instrumente	Physik Instrumente
(PI Shanghai) Co., Ltd.	(PI Shanghai) Co., Ltd.
Shanghai	Beijing
www.pi-china.cn	www.pi-china.cn
SOUTHEAST ASIA	TAIWAN
PI (Physik Instrumente)	Physik Instrumente (PI)
Singapore LLP	Taiwan Ltd.
Singapore	Taipei
www.pi-singapore.sg	www.pi-taiwan.com.tw
For ID / MY / PH / SG /TH / VNM	
KOREA	
PI Korea Ltd. Seoul	



V-308 Voice Coil PIFOC

HIGH DYNAMICS OBJECTIVE FOCUS DRIVE WITH WEIGHT FORCE COMPENSATION

www.pikorea.co.kr

Introducing the V-308 Voice Coil PIFOC

Highly Versatile Z Axis Positioner

From genome sequencing through multiphoton fluorescence microscopy, deep tissue inspection, laser materials processing, and wafer inspection to magnetic tweezers in research ... vertical positioning systems with large travel range, high dynamics, and precision are required in numerous applications in microscopy and industry. With the new V-308 voice coil PIFOC, PI now offers a magnetic direct drive based solution that unites all these requirements.

The stage can be operated with the specially developed C-414 controller (see image) using the proven PIMicroMove user interface. Alternatively, the stage can also be operated with several ACS controllers, namely the SPiiPlus controller supplemented by the NPMpm driver/amplifier for even higher dynamics.



Key Features

Extraordinary Performance for Your Applications

High dynamics

- Accelaration: 8m/s²
- Max. speed: 200 mm/s
- Step-and-settle: <15 ms @100 nm and 250 nm; ±15 nm

Precision in motion

- Min. Step size: 10 nm
- Bidirectional repeatability: 25 nm (on 100 nm travel)
- Position noise: 3 nm



Flexible usage and mounting

- Adapter rings for simple exchange of objective
- Adapter plate for mounting the drive to an optical plate
- Objective holder for horizontal or vertical assembly

Compact design / proven precision components

- Compact outer dimensions: 87 x 77 x 30 mm (H/W/D)
- Crossed roller bearings with anti-creep system
- Reliable voice coil motor
- Proven PIOne linear encoder

User-adjustable weight force compensation: to 1 kg



ACS

4 >> For a complete overview of specifications please click here to visit the product webpage

Key Benefits

High Productivity, Fast, and Reliable Results

EXTENSIVE ACCESSORIES



Adapter rings for objective mounting



Adapter plate for mounting on an optical bench or plate



Holders for vertical or horizontal objective mounting

User-adjustable travel range/stops from 1 mm to 7 mm

2 or 3-Photon Fluorescence Microscopy

even three) excite visible fluorescence in deep tissue layers in the infrared range. By using infrared emissions, a deeper penetration is achieved; this is why a respectively large travel range is required. For example, cell





Deep Tissue Multiphoton Fluorescence Microscopy

The deep tissue multiphoton fluorescence microscopy is a subtype of the 2-photon fluorescence microscopy. Here, Z stacks with a thickness of 6 mm are created – correspondingly large is the travel range required for the Z positioner. This type of microscopy can, for example, be applied in the exploration of the connectome, the interconnection of neurons in the brain. Typically, inverse, confocal microscopes are used with special objectives that have a mass of up to 800 g.

Genome Sequencing

An essential method for genome analysis is the so-called 'squencing-bysynthesis' method. With this method, nucleotide with fluorescent markers are gradually incorporated into a single strand DNA and are so 'synthesized'. If it comes to a recombination of the base pair, a characteristic fluorescent color signal ('occurrence') is generated for each of the four different base pair combinations. These signals are detected microscopically. Since the signals are weak and also very short lived, the detection and imaging process must be very fast. For this purpose, the objective must be precisely aligned at the height of the occurrence along the optical axis, within the depth of field of the lens. Here, the required travel ranges are up to 5 mm.





Digital Slide Scanning

This application is about creating and achieving digital sample images by widely automating the scanning of tissue samples, mainly from routine examinations. These images can then be can be analyzed from experts around the world and samples do not have to be placed repeatedly under the microscope. For example, this is important for following up medical treatments.

The required travel ranges of one to two millimeters are not as large as in other applications; however, high dynamics when focusing is very important.

In two or three photon fluorescence microscopy, two photons (sometimes organs, metabolic processes, etc. were examined.